

CLASSIC SERIES | CORTICAL™

RECOMMENDED STRAIGHT DRILL PROTOCOL

		Drill Diameter [mm]	Ø1.9	Ø2.0	Ø2.8	Ø3.2	Ø3.65
		Drill Speed [RPM]	1200-1500	900-1200	800-1000	500-700	400-700
IMPLANT DIAMETER	Ø4.0	Soft Bone	▼	→	▼		
		Hard Bone	▼	→	▼	→	▼
	Ø5.0	Soft Bone	▼	→	▼		
		Hard Bone	▼	→	▼	→	▼
	Ø6.0	Soft Bone	▼	→	▼		
		Hard Bone	▼	→	▼	→	▼

RECOMMENDED STEP DRILL PROTOCOL

		Drill Diameter [mm]	Ø1.9	Ø2.0	Ø2.8	Ø3.2	Ø3.65
		Drill Speed [RPM]	1200-1500	900-1200	800-1000	500-700	400-700
IMPLANT DIAMETER	Ø4.0	Soft Bone	▼	→	▼		
		Hard Bone	▼	→	▼	→	▼
	Ø5.0	Soft Bone	▼	→	▼		
		Hard Bone	▼	→	▼	→	▼
	Ø6.0	Soft Bone	▼	→	▼		
		Hard Bone	▼	→	▼	→	▼

	Drill to mark osteotomy site		Drill osteotomy to implant		Drill osteotomy partially according to implant
---	------------------------------	---	----------------------------	---	--

The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment. The implants may be loaded for immediate function when good primary stability (above 35 Ncm) has been achieved and with appropriate occlusal loading.