


ONE-PIECE SERIES | MONO™



BONE TYPES	All bone types
DESIGN FEATURES	<ul style="list-style-type: none"> • Tapered thread and tapered core body • Cementable prosthetic portion
CLINICAL BENEFITS	<ul style="list-style-type: none"> • Tissue level implant • Bone condensing • High primary stability • Minimal drilling • Immediate loading

ORDERING INFORMATION



Ø D (mm)	Ø D1 (mm)	Ø D2 (mm)	L (mm)	Ref. No
3.0	1.8	2.0	8	NM-V3008
			10	NM-V3010
			11.5	NM-V3011
			13	NM-V3013
			16	NM-V3016
3.3	1.8	2.1	6	NM-V3306
			8	NM-V3308
			10	NM-V3310
			11.5	NM-V3311
			13	NM-V3313
3.75	1.9	2.5	16	NM-V3316
			6	NM-V3706
			8	NM-V3708
			10	NM-V3710
			11.5	NM-V3711
4.2	1.9	2.8	13	NM-V3713
			16	NM-V3716
			6	NM-V4206
			8	NM-V4208
			10	NM-V4210
5.0	1.9	2.8	11.5	NM-V4211
			13	NM-V4213
			16	NM-V4216
			18	NM-V4218
			6	NM-V5006
			8	NM-V5008
5.0	1.9	2.8	10	NM-V5010
			11.5	NM-V5011
			13	NM-V5013
			16	NM-V5016

ONE-PIECE SERIES | MONO™

RECOMMENDED STRAIGHT DRILL PROTOCOL

Drill Diameter [mm]		Ø1.5	Ø2.0	Ø2.8	Ø3.2	Ø3.65
Drill Speed [RPM]		1500-1200	1200-900	1000-800	700-500	700-400
IMPLANT DIAMETER	Ø3.0	Soft Bone	↓ → 1/3 ↓			
		Hard Bone	↓ → 1/3 ↓			
	Ø3.3	Soft Bone	↓ → 2/3 ↓			
		Hard Bone	↓ → ↓			
	Ø3.75	Soft Bone	↓ → ↓			
		Hard Bone	↓ → ↓ → 1/3 ↓			
	Ø4.2	Soft Bone	↓ → ↓ → 2/3 ↓			
		Hard Bone	↓ → ↓ → 2/3 ↓ → 1/3 ↓			
	Ø5	Soft Bone	↓ → ↓ → 2/3 ↓ → 1/3 ↓			
		Hard Bone	↓ → ↓ → 2/3 ↓ → 2/3 ↓ → 1/3 ↓			








RECOMMENDED STEP DRILL PROTOCOL

Drill Diameter [mm]		Ø1.9	Ø2.0	Ø2.8	Ø3.2	Ø3.65
Drill Speed [RPM]		1500-1200	1200-900	1000-800	700-500	700-400
IMPLANT DIAMETER	Ø3.0	Soft Bone	↓			
		Hard Bone	↓ → 1/3 ↓			
	Ø3.3	Soft Bone	↓ → 1/3 ↓			
		Hard Bone	↓ → 2/3 ↓			
	Ø3.75	Soft Bone	↓ → 2/3 ↓			
		Hard Bone	↓ → ↓ → 1/3 ↓			
	Ø4.2	Soft Bone	↓ → ↓ → 2/3 ↓			
		Hard Bone	↓ → ↓ → 2/3 ↓ → 2/3 ↓			
	Ø5	Soft Bone	↓ → ↓ → 2/3 ↓ → 2/3 ↓			
		Hard Bone	↓ → ↓ → 2/3 ↓ → 2/3 ↓ → 2/3 ↓			

The recommended drill protocol procedure should not replace the dentist's/surgeon's judgment. The implants may be loaded for immediate function when good primary stability (above 35 Ncm) has been achieved and with appropriate occlusal loading.

 Drill to mark osteotomy site	 Drill osteotomy to implant	 Drill osteotomy partially according to implant
--	--	--

COMPONENTS

						
NMCD2315 Drill 1.5mm Ø	NM-X1620 Motor Driver 2.15mm Ø	NM-X1018 Ratchet Driver 2.15mm Ø	NM-X1019 Ratchet Driver 2.15mm Ø	NM-X1720 Ratchet Driver 2.15mm Ø	NM-T6202 Analog	NM-T4420 Transfer