


CEMENT RETAINED RESTORATIONS | ANGLED ABUTMENTS



DESIGN FEATURES	<ul style="list-style-type: none"> • Angled abutments: 15° and 25° • Wide range of diameters and lengths • Prefabricated with various collar heights
CLINICAL BENEFITS	<ul style="list-style-type: none"> • Cement retained abutments • For dentist/technician customization • Wide range of emergence profile abutment options
MATERIAL	Titanium Alloy (Ti-6Al-4V ELI)

ORDERING INFORMATION



ANATOMIC				
Angle 15°	Angle 25°	Ø D (mm)	L (mm)	L1 (mm)
NM-A2401	NM-A2601	5.4	9.5	1
NM-A2402	NM-A2602	5.4	10.5	2
NM-A2403	NM-A2603	5.4	11.5	3
NM-A2404	NM-A2604	5.4	12.5	4

STANDARD			
Angle 15°	Angle 25°	Ø D (mm)	L (mm)
NM-A3209	NM-A3409	4.5	9
NM-A3211	NM-A3411	4.5	11

NARROW TOP			
Angle 15°	Angle 25°	Ø D (mm)	L (mm)
NM-A2809	NM-A3009	4.5	9

NARROW			
Angle 15°	Angle 25°	Ø D (mm)	L (mm)
NM-A3609	NM-A3809	4	9



NM-S8324
Included with all abutments & available separately