


# CEMENT RETAINED RESTORATIONS | TEMPORARY ABUTMENTS



DESIGN FEATURES	<ul style="list-style-type: none"> <li>• Sterilizable</li> <li>• Connects to the implant with a retentive screw</li> <li>• High compressive strength</li> </ul>
CLINICAL BENEFITS	<ul style="list-style-type: none"> <li>• Acrylic retained abutment</li> <li>• Does not create artifacts in X-ray irradiation procedures, magnetic resonance imaging (MRI), and computed tomography (CT)</li> <li>• Does not contain metal additives to prevent ion exchange in the mouth</li> <li>• Easy and fast customization by the dentist or technician</li> </ul>
MATERIAL	Peek

## ORDERING INFORMATION



Straight	L (mm)	Angle 15°	Angle 25°	L (mm)	L1 (mm)
	9			9.5	1
NM-C6001		NM-C5001	NM-C5201		
	10			10.5	2
NM-C6002		NM-C5002	NM-C5202		
	11			11.5	3
NM-C6003		NM-C5003	NM-C5203		



NM-S8324  
Included with all abutments  
& available separately