


ACCESSORIES | CLIP-ON TRANSFER




DESIGN FEATURES	<ul style="list-style-type: none"> • Retentive split pin • Clips in to the implant without the need of a screw
CLINICAL BENEFITS	<ul style="list-style-type: none"> • Enables closed tray impression transfer • Especially useful in limited opening situations • Shortens working time. No need to screw in or screw out
MATERIAL	Stainless Steel

ORDERING INFORMATION



Transfer	Ø D (mm)	L (mm)
	4.7	9
NM-T3409		
	4.7	13
NM-T3413		

ACCESSORIES | IMPLANT ANALOG



Analog	Ø D (mm)	L (mm)
	4	12.7
NM-T6004		
	5	12.7
NM-T6005		

