

Pterygoid Surgical Set

READY...SET...PTERY-GO!

Introducing the Noris **Pterygoid Surgical kit**.

A one of a kind dedicated kit for the surgical placement of Noris pterygoid implants.

The kit has osteotomes and surgical drills to establish the path for the pterycore or pteryfit implants.

The surgical screwdriver, torque ratchet, or motor mount drivers enable easy and precise placement of the implant in to place.


Finally, the manual or motor hex drivers facilitate the attachment of the appropriate prosthetic components onto the Noris pterygoid implant.

Get set for predictable full arch maxillary surgery with the Noris Pterygoid kit.



PTERYGOID SURGICAL SET

NM-X2250



LBL061 Rev.2022-1

SURGICAL SET FOR PTERYGOID NM-X2250 CONTAINS:

PR-X2250	Pterygoid Surgical Box	NMHX1015	Motor Mount Hex 2.42, L28.0mm, Self Holding, SS.
NMCD2320	Drill, \varnothing 2.0, L25mm, DLC, SS.	MN-X1215	Driver Hex 1.25mm, Length 15.0mm, SS.
NMCD4032	Conical Drill, \varnothing 2.0- \varnothing 3.2mm, L25mm, DLC, SS.	NM-X7011	Hand driver Star-Hex 1.25, Length 20.0mm, SS.
NM-D7432	Drill, \varnothing 3.2mm, Length 43.0mm, SS.	NM-X1010	Motor Mount Hex 1.25, Length 35.0mm, SS.
NM-X3001	Osteotom 2.0 mm, SS	NM-D3412	Drill Extension, SS.
NM-X3002	Osteotom 3.0 mm, SS	NM-X1021	Torque ratchet
NMHX2615	Driver Hex 2.42mm, L15.0mm, Self Holding, SS	NM-X1023	Surgical screw driver long
NMHX2640	Driver Hex 2.42mm, L40.0mm, Self Holding, SS		

