

KITS | KIT QUIRÚRGICO PREMIUM

NMCX2111 / NMCX2112



EL KIT NMCX2111 CONTIENE:

Fresa, \varnothing 2.0mm	NMCD1220
Fresa, \varnothing 2.8mm	NMCD1228
Fresa, \varnothing 3.2mm	NMCD1232
Fresa, \varnothing 3.65mm	NMCD1236
Fresa, \varnothing 4.2mm	NMCD1242
Fresa, \varnothing 5.2mm	NMCD1252

EL KIT NMCX2112 (FRESAS CÓNICAS) CONTIENE:

Fresa, \varnothing 1.5mm	NMCD2315
Fresa, \varnothing 1.8-2.4mm	NMCD3018
Fresa, \varnothing 2.0-3.2mm	NMCD3020
Fresa, \varnothing 2.5-3.7mm	NMCD3025
Fresa, \varnothing 2.7-4.5mm	NMCD3028
Fresa, \varnothing 3.1-5.5mm	NMCD3031

AMBOS KITS CONTIENEN:

Fresa avellanadora, \varnothing 3.8-4.2mm	NMCD1034	Destornillador de motor \varnothing 1.25, L 20mm	NM-X1008
Fresa avellanadora, \varnothing 5.0-6.0mm	NMCD1056	Insertor de motor \varnothing 2.42, L 20mm	NMHX1014
Fresa marcadora	NM-D3410	Insertor de motor \varnothing 2.42, L 28mm	NMHX1015
Prolongador de Fresas	NM-D3412	Carraca dinamométrica	NM-X1021
Destornillador de carraca \varnothing 1.25mm, L 7mm	NM-X1207	Paralelizador quirúrgico recto, L 10mm - 2 Pcs.	NM-X1026
Destornillador de carraca \varnothing 1.25mm, L 15mm	NM-X1215	Paralelizador quirúrgico recto, L 16mm - 2 Pcs.	NM-X1027
Insertor manual \varnothing 2.42mm, L 7mm	NMHX2607	Medidor de profundidad	NM-X1025
Insertor manual \varnothing 2.42mm, L 15mm	NMHX2615	Destornillador manual \varnothing 1.25, L 14mm	NM-X1007