


CLASSIC SERIES | TUFF PRO™



BONE TYPES	All bone types
PROSTHETICS PLATFORM	Internal hex
DESIGN FEATURES	<ul style="list-style-type: none"> • Condensing variable threads design • Apically tapered threads and tapered core body • Double threads with large step • Machined surface coronal portion • Double flutes
CLINICAL BENEFITS	<ul style="list-style-type: none"> • Self tapping • High primary stability • Minimal drilling • Fast insertion – optimal for soft bone • Immediate loading - suitable for extraction sites

ORDERING INFORMATION



Ø D (mm)	Ø D0 (mm)	Ø D1 (mm)	L (mm)	Ref. No
3.3	1.5	2.6	8	NMAF3308
			10	NMAF3310
			11.5	NMAF3311
			13	NMAF3313
			16	NMAF3316
3.75	1.8	3.1	8	NMAF3708
			10	NMAF3710
			11.5	NMAF3711
			13	NMAF3713
			16	NMAF3716
4.2	2.1	3.5	6	NMAF4206
			8	NMAF4208
			10	NMAF4210
			11.5	NMAF4211
			13	NMAF4213
5.0	2.7	4.5	16	NMAF4216
			6	NMAF5006
			8	NMAF5008
			10	NMAF5010
			11.5	NMAF5011
			13	NMAF5013
			16	NMAF5016

Cover Screw Included in all Internal Hex implants



NM-S5023

