

KITS | KIT QUIRÚRGICO ZYGOMATIC

NM-X2118



EL KIT QUIRÚRGICO ZYGOMATIC NM-X2118 CONTIENE:

Fresa {1} $\varnothing 2.0/\varnothing 2.5/\varnothing 3.0$, L60mm	NM-D7020
Fresa {2} $\varnothing 2.8/\varnothing 3.0$, L60mm	NM-D7028
Fresa {3} $\varnothing 2.8/\varnothing 3.2$, L60mm	NM-D7032
Fresa {1} $\varnothing 2.0/\varnothing 2.5/\varnothing 3.0$, L80mm	NM-D7120
Fresa {2} $\varnothing 2.8/\varnothing 3.0$, L80mm	NM-D7128
Fresa {3} $\varnothing 2.8/\varnothing 3.2$, L80mm	NM-D7132
Fresa, $\varnothing 4.2$ mm, L30mm, grano fino	NM-D7201
Fresa, $\varnothing 4.2$ mm, L30mm, grano medio	NM-D7202
Fresa, $\varnothing 4.2$ mm, L30mm, grano grueso	NM-D7203
Medidor de profundidad de 35mm a 60mm	NM-X1028
Insertor manual hex. 2.4mm, L20mm	NMHX2620
Destornillador de carraca para prótesis 1.25mm, L10mm	NM-X1210
Destornillador manual hex. de estrella 1.25, L14mm	NM-X7007